



# Animal Hyperbaric

## Animal Hyperbaric Chamber range

### The best animal system in hyperbaric technology!

✓ 100% Natural Zero side effects ✓ Improved oxygen levels ✓ Acute & Chronic health conditions ✓ Faster Recovery

## Vetsmart Clear acrylic Hyperbaric Chamber

Have the best possible treatment available for your pet today!

In the animal hyperbaric chamber, pets lie on a soft blanket and rest or sleep while the oxygen goes to work on wounds, swelling, burns and other injuries or illnesses. The pets are comfortable and relaxed during dog/ cat Hyperbaric Oxygen Therapy - an important element to faster recovery. The total veterinary HBOT treatment time is from 1 to 2 hours, and is usually repeated twice a day. Treatments continue until the doctors see a remarkable improvement. When you pet is successfully starting to use the affected limb, or is gaining strength and function, the animal Hyperbaric oxygen chamber treatments are discontinued. All systems options include Air cooler, Air compressor, High pressure medical oxygen concentrator. External pressure gages. Manual pressure relief valves. Emergency pressure relief valves.



High Pressure  
medical oxygen  
concentrator  
specifications:

Voltage	110/220V
Oxygen Flow	1-10L High Pressured
Power	480W
Weight	30kg

Air Compressor  
Specifications:

Voltage	110/220V
Air Flow	72L/Minute
Power	480W
Weight	18kg

Air Cooler  
Specifications:

Voltage	110/220V
Capability	150L/Minute
Pressure	0.8mPa
Power	220W
Weight	5kg

Chamber Specifications:

Size	80cm width x 120cm long
Pressure	1,3 ATA
Weight	
Material	Clear acrylic
Monitoring	Temperature
	Humidity
	Pressure
	Oxygen Flow