



Recovery Technology Physiotherapy Technology

The best professional choice for private or commercial use

✓ Muscle Pathology ✓ Bone Injuries ✓ Inflammation and nerve damage



3 in 1 Smart TECAR wave

Used mostly in physiotherapy healing, Smart Tecar Wave is a 3 in 1 therapy machine desktop device, combining Tecar therapy with shockwave therapy and Electro muscle stimulation therapy to speed up the healing process. First, electrical stimulation can be used to relax and stimulate muscle and surface muscles, and then shock wave therapy can be used to hit tendons and deep muscles at a deeper level, and Tecar therapy with 448Khz and 300W can be used to repair the treatment site. You can adjust according to the actual treatment.

Specifications:

Power	1000 watt and 220 Volts
Technology:	Tecar Electro muscle stimulation and Shockwave
Tecar Output power	300watt
Tecar Frequency	300-448KHz
Electro Stimulation Intensity	1-20

Electro Stimulation frequency	1-30Hz
Electromagnetic Shockwave Therapy	5-200Mj(200Mj=4Bar)
Frequency	1-16Hz
Treatment heads	7PCS (2 specific for ED treatment)

Head sizes	20/40/60 mm
Screen	10,4 Inch touch screen
Size	52cm x 51cm x 46cm
Gross weight	30 kg
Electro Stimulation Intensity	1-20



RET Tips



CET Tips



Return Pad



Self vacuum cup



1# ~ 7# heads (number 6 and 7 for ED therapy)