



Aesthetic Beauty Rejuvenation Therapy

The best aesthetic beauty solutions right here!

✓ Pelvic floor therapy ✓ Electro Muscle Stimulation ✓ Rejuvenation Therapy

Pelvic Floor stimulator therapy Chair

6 mode electro muscle stimulation system for both men and woman



System
Specifications:

| | |
|------------|----------------------------------|
| Voltage | 110/220V |
| Dimensions | Length 49cm x Width 48cm |
| Frequency | 50-60 HZ |
| Effect | 12000 Contractions in 30 Minutes |
| Material | ABS and Silica Gel |
| Weight | 5kg |

Treatment period: In this phase , the pelvic floor muscles are treated with intermediate frequency current and thermal energy. It is divided into three phases according to the program mode . For treatment , according to the individual's endurance , the muscle relaxation and contraction time of each stage is different .

The propulsion period of the model is 7 days , and the intensity should be adjusted appropriately according to the individual' s endurance .

Stages for each phase of treatment includes Initial stage, Mid-term, Late stage, Maintenance period and Nursing phase

Treatment is successful for Urinary incontinence, Anal disfunctions, Vaginal wall treatment, Prostate therapy, Impotence therapy and rectification of premature Ejaculation therapy.

Pelvic chair helps with active regeneration of collagen in pelvic floor tissues, strengthens pelvic floor muscles, Improve privacy issues, such as urinary incontinence, penetration. Achieve vagina tightening.

The total treatment time in one day may not exceed 1.5hours, one individual body part does not exceed 1 hour treatment time at once. Optimal results 6 sessions

Non-invasive treatment with no radiation, pain or side effects and no labour for your therapist

